

Introduction

Large-scale expansion of human pluripotent stem cells is an evolving field with promising applications in cell therapy and regenerative medicine. Despite significant advances, economic viability of large-scale stem cell culture bioprocesses remains challenging, hampering translation. Cost of cell **culture media** is a major barrier.

The current work sought to answer the following question: **How does optimizing culture medium composition, perfusion culture, and online monitoring of glucose & lactate concentrations help to reduce culture media associated costs?** To answer, scientists at IST performed experiments using the DASbox[®] Mini Bioreactor System, culturing iPSCs as free-floating scaffold-free aggregates. Central to this work is HiDef[®] S8, a new chemically defined, animal-free suspension medium developed specifically for the efficient expansion and maintenance of human pluripotent stem cells. Daily samples were drawn from the bioreactors for cell counting, viability measurements, and immunocytochemistry was performed for flow cytometry.

Background

Stirred-tank bioreactors are ideal culture platforms for the large-scale culture of **stem cells**, because they:

- Integrate monitoring and automation
- Better mimic a 3D cellular environment
- Facilitate high cell densities
- Are highly scalable

Human induced pluripotent stem cells (iPSCs) are a promising cell type for advancing clinical applications.


This work tested four culture media for iPSC expansion: **mTeSR1**, **E8** and two formulations from **Defined Bioscience**: **HiDef-B8** and **HiDef-S8**.

Materials and Methods


Experiments were performed using Eppendorf's DASbox[®] Mini Bioreactor System, culturing iPSCs as freefloating scaffold-free 3D aggregates:

- Inoculation density of 100,000 iPSCs / mL
- Working volume of 80 mL
- Media exchange via continuous perfusion
- Monitoring of pH, DO, Temp, & glucose & lactate


Cell counts, viability staining for Calcein AM, and immunocytochemistry for FC performed daily.





Results



mTeSR1 and **B8** led to fold expansions of 33× in 7 days, whereas **E8** led to a fold expansion of 18× in 7 days and **HiDef-S8** led to a fold expansion of 39× in 6 days. **mTeSR1** and **HiDef-S8** presented improved cell survival (38% and 37%, respectively). **HiDef-S8**'s growth rate was significantly faster than when using other culture media.




HiDef-S8 expanded aggregates present a hollow morphology (above). All aggregates stain for Calcein AM throughout expansion, indicating viability (below).



For pluripotency maintenance, only **mTeSR1** and **HiDef-S8** preserved high expression of key pluripotency markers (>95%). **E8** and **B8** both led to a significant decrease in TRA-1-60 expression at the end of culture.

FC Marker	mTeSR1	E8	B8	HiDef-S8
OCT4	97.0	96.5	86.2	95.2
SOX2	98.0	99.0	96.0	99.3
TRA160	97.8	72.9	46.4	97.4
SSEA-4	98.2	97.9	97.8	98.4

mTeSR1 and **HiDef-S8** led to greater iPSC (TRA-1-60+) yields than **E8** and **B8**. Despite the similar performance of **mTeSR1** and **HiDef-S8** in this metric, **HiDef-S8**'s significantly lower cost and faster growth rate led to >34% cost reduction.



iPSC expanded in **HiDef-S8** over multiple iterations, and multiple passages maintain aggregate diameter below 400 μm, present normal karyotype, and readily differentiated into all three germ layers in adherent culture post-3D expansion.

Conclusions

For all media conditions, aggregates presented a hollow morphology and stained positive for Calcein AM viability marker throughout 6-7 days. **HiDef-S8** led to a fold expansion of 39× in 6 days. **HiDef-B8** and **TeSR** both led to fold expansions of 33× in 7 days, and **E8** led to a fold expansion of only 18× in 7 days. **TeSR** and **S8** presented improved cell survival.

S8 supported faster expansion and consumed volumes less medium. **HiDef-S8** and **TeSR** led to much greater TRA-1-60+ yields over 6 and 7 days, respectively, than either **E8** or **B8**. Despite the similar performance of **TeSR** and **S8** in this metric, **S8**'s **lower cost per volume** and **faster stable growth rate** resulted in **substantially lower medium cost**. In terms of pluripotency maintenance, only **HiDef-S8** and **TeSR** preserved high expression of key pluripotency markers (>95%) after 6 days and 7 days, respectively. Notably, for extracellular vesicle production it is important that the culture medium contain minimal particles. Nanoparticle analysis showed that **TeSR contained substantially higher concentrations of particles** across multiple lots compared to **HiDef-S8**, **HiDef-B8**. **HiDef-S8**, a new suspension medium specifically optimized to support high-density PSC expansion, led to a 34% reduction in culture media associated costs comparatively to **TeSR**. Perfusion and glucose/lactate monitoring also decreases costs, by facilitating an optimal medium exchange rate.

Future Work

Future work will focus on manufacturing cGMP-compliant **HiDef-S8** to support clinical applications. Defined Bioscience aims to expand their 3D culture platform to enable differentiation of human PSCs in suspension to further enhance scalability and efficiency for large-scale cell production.

Funding

- Fundação para a Ciência e a Tecnologia;
- I.P. granted to CardioWheel (PTDC/EQU-EQU/29653/2017) and SMART (PTDC/EQU-EQU/3853/2020)
- I.P. granted to iBB (UIDB/04565/2020 and UIDP/04565/2020)

Acknowledgments