

Transforming Growth Factor Beta-3 in Suspension Culture Media for Induced Pluripotent Stem Cells

Jannah Jitkaroon^{1,2,3}, Gerardo Aguirre Castillo¹, Rebecca Pierce¹, Maria Mozo¹, Jerome V. Karpiak¹, Steven D. Rees^{1*}

¹Defined Bioscience Inc, San Diego CA, 6404 Nancy Ridge Dr., San Diego, CA, 92121 *steve@definedbioscience.com

²California Institute for Regenerative Medicine

³California State University San Marcos




Introduction


Growth factors used in stem cell culture media account for approximately 80% of cost, making it one of the highest cost drivers of cell therapy and in-vitro biomedical research. Transforming growth factor beta (TGF β) is a family of proteins that promote cell growth, embryogenesis, and differentiation and is essential in pluripotent stem cell medium formulations to support survival and proliferation. TGF β -1 is most used in medium formulations, however, its signal promotes spontaneous differentiation and culture heterogeneity over extended expansion. TGF β -3 activates the same SMAD2/3 pathway but supports pluripotency maintenance at a twenty-fold lower working concentration, resulting in long-term culture robustness and reduced costs. We expressed recombinant TGF β -3 in *E. coli* and purified it using chromatography. Then, we compared effects in suspension culture by analyzing aggregate morphology, size, expansion, and pluripotency.

Methods

TGFB3 Production



Suspension Culture



Results




Figure 1: Images of iPSCs on day 4 of suspension culture (4X magnification). Similar morphology of aggregates (round, defined borders, no hollowness) in S8 (A), S8 with lower TGF β -3 (B), and S8 with in-house TGF β -3 (C) was seen.




Figure 2: Aggregate Size and Expansion Analysis. iPSCs were cultured for four days in Defined Bioscience's suspension culture medium (HiDef® S8®). The width of individual aggregates were measured using ImageJ software on Day 4 (A). The fold change was also assessed (B). Decreasing TGF β -3 yielded significantly larger aggregates (p<0.05) and highest expansion.




Figure 3: Pluripotency Analysis. Day 4 iPSC aggregates were stained with SSEA-4 conjugated with APC and TRA-1-60 conjugated with PE antibodies. Samples were analyzed with the Novocyte flow cytometer. Results indicate similar levels of SSEA-4 and TRA-1-60 expression across S8 (A), S8 with lower TGF β -3 (B), and S8 with in-house TGF β -3 (C).

Conclusions and Future Directions

Conclusions:

- TGF β -3 was produced via recombinant expression in *E. Coli*
- Aggregates across all conditions had similar morphology
- Lowering TGF β -3 maintains similar levels of pluripotency expression while having the highest fold change
- Using in-house TGF β -3 did not have an impact on aggregate size, while maintaining pluripotency

Future Directions:

- Assess iPSC aggregate differentiation potential
- Optimization of suspension culture by modulating shear stress
- Optimizing production and purification of in-house TGF β -3
- Serial expansion of iPSC aggregates
- Frozen down cell aggregates for easy thaw directly into suspension

References

Kuo, H.-H., Gao, X., DeKeyser, J.-M., Fetterman, K. A., Pinheiro, E. A., Weddle, C. J., Fonoudi, H., Orman, M. V., Romero-Tejeda, M., Jouni, M., Blancard, M., Magdy, T., Epting, C. L., George, A. L., & Burrige, P. W. (2020). *Negligible-cost and weekend-free chemically defined human iPSC culture. Stem Cell Reports*, 14(2), 256–270. <https://doi.org/10.1016/j.stemcr.2019.12.007>

Acknowledgments

